

Product datasheet for **TB409649**

HDAC6 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HDAC6 (NM_006044) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209649
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	HDAC6 Mouse Monoclonal Antibody [Clone ID: OTI4C5] (TA502030)
ACCN:	NM_006044 , NP_006035
Synonyms:	CPBHM; HD6; JM21; PPP1R90
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_006035
Locus ID:	10013
Cytogenetics:	Xp11.23
Protein Families:	Druggable Genome, Transcription Factors



[View online »](#)