

Product datasheet for **TB409729**

Hormone sensitive lipase (LIPE) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of LIPE (NM_005357) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209729
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_005357 , NP_005348
Synonyms:	AOMS4; FPLD6; HSL; LHS; REH
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005348
Locus ID:	3991
Cytogenetics:	19q13.2
Protein Pathways:	Insulin signaling pathway



[View online »](#)