

Product datasheet for **TB409772**

RPL13 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RPL13 (NM_000977), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209772
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_000977</u> , <u>NP_000968</u>
Synonyms:	BBC1; D16S44E; D16S444E; L13; SEMDIST
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_000968</u>
Locus ID:	6137
Cytogenetics:	17p11.2
Protein Families:	Druggable Genome
Protein Pathways:	Ribosome



[View online »](#)