

Product datasheet for **TB409776**

RNA Polymerase II p14.5 (POLR2I) CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of POLR2I (NM_006233) in HEK293T cells paraffin embedded controls for ICC/IHC staining |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC209776 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100) |
| ACCN: | NM_006233 , NP_006224 |
| Synonyms: | hRPB14.5; RPB9 |
| Storage: | Room Temperature, or 2-8°C for long term storage |
| Stability: | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_006224 |
| Locus ID: | 5438 |
| Cytogenetics: | 19q13.12 |
| Protein Families: | Transcription Factors |
| Protein Pathways: | Huntington's disease, Metabolic pathways, Purine metabolism, Pyrimidine metabolism, RNA polymerase |



[View online »](#)