

Product datasheet for **TB409872**

GPR126 (ADGRG6) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ADGRG6 (NM_001032395), transcript variant b2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC209872
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_001032395 , NP_001027567
Synonyms:	APG1; DREG; GPR126; LCCS9; PR126; PS1TP2; STQTL1; VIGR
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001027567
Locus ID:	57211
Cytogenetics:	6q24.2
Protein Families:	Druggable Genome, GPCR, Transmembrane



[View online »](#)