

Product datasheet for **TB409987**

RPL10A CytoSection

Product data:

| | |
|---------------------------------------|---|
| Product Type: | CytoSections |
| Description: | Transient overexpression of RPL10A (NM_007104) in HEK293T cells paraffin embedded controls for ICC/IHC staining |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC209987 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100) |
| Target Detection Antibodies: | RPL10A Mouse Monoclonal Antibody [Clone ID: OTI2G9] (TA809519) |
| ACCN: | NM_007104 , NP_009035 |
| Synonyms: | Csa-19; CSA19; L10A; NEDD6 |
| Storage: | Room Temperature, or 2-8°C for long term storage |
| Stability: | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_009035 |
| Locus ID: | 4736 |
| Cytogenetics: | 6p21.31 |
| Protein Pathways: | Ribosome |



[View online »](#)