

Product datasheet for **TB410001**

S100 Calcium Binding Protein A13 (S100A13) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of S100A13 (NM_001024210), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210001
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_001024210</u> , <u>NP_001019381</u>
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001019381</u>
Locus ID:	6284
Cytogenetics:	1q21.3
Protein Families:	Druggable Genome



[View online »](#)