

## Product datasheet for **TB410007**

### XCL2 CytoSection

#### Product data:

|                                       |                                                                                                                                                                                                                         |
|---------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:                         | CytoSections                                                                                                                                                                                                            |
| Description:                          | Transient overexpression of XCL2 (NM_003175) in HEK293T cells paraffin embedded controls for ICC/IHC staining                                                                                                           |
| Species:                              | Human                                                                                                                                                                                                                   |
| Expression Host:                      | HEK293T                                                                                                                                                                                                                 |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC210007                                                                                                                                                                                                  |
| Tag:                                  | C-MYC/DDK                                                                                                                                                                                                               |
| Detection Antibodies:                 | Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)                                                                                                                                                         |
| ACCN:                                 | <u><a href="#">NM_003175</a></u> , <u><a href="#">NP_003166</a></u>                                                                                                                                                     |
| Synonyms:                             | SCM-1b; SCM1B; SCYC2                                                                                                                                                                                                    |
| Storage:                              | Room Temperature, or 2-8°C for long term storage                                                                                                                                                                        |
| Stability:                            | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.                                                                                                                 |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.                                                                                                                   |
| RefSeq:                               | <u><a href="#">NP_003166</a></u>                                                                                                                                                                                        |
| Locus ID:                             | 6846                                                                                                                                                                                                                    |
| Cytogenetics:                         | 1q24.2                                                                                                                                                                                                                  |
| Protein Families:                     | Druggable Genome, Secreted Protein                                                                                                                                                                                      |
| Protein Pathways:                     | Chemokine signaling pathway, Cytokine-cytokine receptor interaction                                                                                                                                                     |



[View online »](#)