

Product datasheet for **TB410199**

MMP16 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MMP16 (NM_005941) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210199
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_005941 , NP_005932
Synonyms:	C8orf57; MMP-X2; MT-MMP2; MT-MMP3; MT3-MMP
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005932
Locus ID:	4325
Cytogenetics:	8q21.3
Protein Families:	Druggable Genome, Protease, Secreted Protein, Transmembrane



[View online »](#)