

## Product datasheet for **TB410328**

### UBE2E3 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of UBE2E3 (NM_006357), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210328
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	UBE2E3 Mouse Monoclonal Antibody [Clone ID: OTI7E8] (TA800082)
ACCN:	<u><a href="#">NM_006357</a></u> , <u><a href="#">NP_006348</a></u>
Synonyms:	UBCH9; UbcM2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_006348</a></u>
Locus ID:	10477
Cytogenetics:	2q31.3
Protein Pathways:	Ubiquitin mediated proteolysis



[View online »](#)