

## Product datasheet for **TB410370**

### PAX4 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PAX4 (NM_006193) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210370
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	PAX4 Mouse Monoclonal Antibody [Clone ID: OTI2F3] (TA801976)
ACCN:	<a href="#">NM_006193</a> , <a href="#">NP_006184</a>
Synonyms:	KPD; MODY9
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_006184</a>
Locus ID:	5078
Cytogenetics:	7q32.1
Protein Families:	Embryonic stem cells, ES Cell Differentiation/IPS, Transcription Factors
Protein Pathways:	Maturity onset diabetes of the young



[View online »](#)