

Product datasheet for **TB410475**

SELS (SELENOS) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SELENOS (NM_203472), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210475
Tag:	Tag Free
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_203472</u> , <u>NP_982298</u>
Synonyms:	AD-015; ADO15; SBB18; SELS; SEPS1; VIMP
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_982298</u>
Locus ID:	55829
Cytogenetics:	15q26.3
Protein Families:	Druggable Genome



[View online »](#)