

## Product datasheet for **TB410649**

### Ku80 (XRCC5) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of XRCC5 (NM_021141) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210649
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u>NM_021141</u> , <u>NP_066964</u>
Synonyms:	KARP-1; KARP1; KU80; Ku86; KUB2; NFIV
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_066964</u>
Locus ID:	7520
Cytogenetics:	2q35
Protein Families:	Druggable Genome, Stem cell - Pluripotency
Protein Pathways:	Non-homologous end-joining



[View online »](#)