

Product datasheet for **TB410707**

Visfatin (NAMPT) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NAMPT (NM_005746) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210707
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_005746 , NP_005737
Synonyms:	1110035O14Rik; PBEF; PBEF1; VF; VISFATIN
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_005737
Locus ID:	10135
Cytogenetics:	7q22.3
Protein Families:	Druggable Genome
Protein Pathways:	Nicotinate and nicotinamide metabolism



[View online »](#)