

Product datasheet for **TB410753**

Coronin 1a (CORO1A) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CORO1A (NM_007074), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210753
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Coronin 1a (CORO1A) Mouse Monoclonal Antibody [Clone ID: OT11A5] (TA501429)
ACCN:	<u>NM_007074</u> , <u>NP_009005</u>
Synonyms:	CLABP; CLIPINA; HCORO1; IMD8; p57; TACO
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_009005</u>
Locus ID:	11151
Cytogenetics:	16p11.2



[View online »](#)