

Product datasheet for **TB410815**

Retinoic Acid Receptor gamma (RARG) CytoSection

Product data:

| | |
|---------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CytoSections |
| Description: | Transient overexpression of RARG (NM_000966), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC210815 |
| Tag: | C-MYC/DDK |
| Detection Antibodies: | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569) |
| Target Detection Antibodies: | Retinoic Acid Receptor gamma (RARG) Mouse Monoclonal Antibody [Clone ID: OTI3G1] (TA810381) |
| ACCN: | NM_000966 , NP_000957 |
| Synonyms: | NR1B3; RARC |
| Storage: | Room Temperature, or 2-8°C for long term storage |
| Stability: | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed. |
| Preparation: | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note: | This product is for research use only and is not approved for use in humans or in clinical diagnosis. |
| RefSeq: | NP_000957 |
| Locus ID: | 5916 |
| Cytogenetics: | 12q13.13 |
| Protein Families: | Druggable Genome, Nuclear Hormone Receptor, Transcription Factors |



[View online »](#)