

## Product datasheet for **TB410897**

### Macrophage Inflammatory Protein 3 (CCL23) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CCL23 (NM_145898), transcript variant CKbeta8, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210897
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<a href="#">NM_145898</a> , <a href="#">NP_665905</a>
Synonyms:	CK-BETA-8; Ckb-8; Ckb-8-1; CKb8; hmrb-2a; MIP-3; MIP3; MPIF-1; SCYA23
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_665905</a>
Locus ID:	6368
Cytogenetics:	17q12
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Chemokine signaling pathway, Cytokine-cytokine receptor interaction



[View online »](#)