

Product datasheet for **TB410931**

Axin 2 (AXIN2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of AXIN2 (NM_004655) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC210931
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_004655 , NP_004646
Synonyms:	AXIL; ODCRCS
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_004646
Locus ID:	8313
Cytogenetics:	17q24.1
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Induced pluripotent stem cells
Protein Pathways:	Basal cell carcinoma, Colorectal cancer, Endometrial cancer, Pathways in cancer, Wnt signaling pathway



[View online »](#)