

Product datasheet for **TB411299**

Casein Kinase 2 beta (CSNK2B) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CSNK2B (NM_001320) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC211299
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u>NM_001320</u> , <u>NP_001311</u>
Synonyms:	CK2B; CK2N; Ckb1; Ckb2; CSK2B; G5A; POBINDS
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001311</u>
Locus ID:	1460
Cytogenetics:	6p21.33
Protein Families:	Druggable Genome
Protein Pathways:	Adherens junction, Tight junction, Wnt signaling pathway



[View online »](#)