

Product datasheet for **TB411541**

RHOBTB1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RHOBTB1 (NM_198225), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC211541
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_198225</u> , <u>NP_937868</u>
Synonyms:	1700008H16Rik; 3110048G13Rik; AI173445; AI573858; AV350930; KIAA0740, MGC33059, MGC33841; MGC107538; Rho-related BTB domain containing 1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_937868</u>
Locus ID:	9886
Cytogenetics:	10q21.2



[View online »](#)