

## Product datasheet for **TB411690**

### Kv beta 2 (KCNAB2) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of KCNAB2 (NM_003636), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC211690
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_003636</u> , <u>NP_003627</u>
Synonyms:	AKR6A5; HKvbeta2; HKvbeta2.1; HKvbeta2.2; KCNA2B; KV-BETA-2
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_003627</u>
Locus ID:	8514
Cytogenetics:	1p36.31
Protein Families:	Druggable Genome, Ion Channels: Other



[View online »](#)