

Product datasheet for **TB412361**

TAF4B CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TAF4B (NM_005640) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC212361
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u>NM_005640</u> , <u>NP_005631</u>
Synonyms:	SPGF13; TAF2C2; TAFII105
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_005631</u>
Locus ID:	6875
Cytogenetics:	18q11.2
Protein Pathways:	Basal transcription factors, Huntington's disease



[View online »](#)