

## Product datasheet for **TB412460**

### IL1RA (IL1RN) CytoSection

#### Product data:

|                                       |   |
|---------------------------------------|---|
| Product Type:                         | CytoSections  |
| Description:                          | Transient overexpression of IL1RN (NM_000577), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining  |
| Species:                              | Human   |
| Expression Host:                      | HEK293T   |
| Expression cDNA Clone or AA Sequence: | TrueORF Clone RC212460  |
| Tag:                                  | C-MYC/DDK   |
| Detection Antibodies:                 | DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)  |
| Target Detection Antibodies:          | IL1RA (IL1RN) Mouse Monoclonal Antibody [Clone ID: OTI3B1] (TA803396)   |
| ACCN:                                 | <a href="#">NM_000577</a> , <a href="#">NP_000568</a>   |
| Synonyms:                             | DIRA; ICIL-1RA; IL-1ra; IL-1ra3; IL-1RN; IL1F3; IL1RA; IRAP; MVCD4  |
| Storage:                              | Room Temperature, or 2-8°C for long term storage  |
| Stability:                            | Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.   |
| Preparation:                          | HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin. |
| Note:                                 | This product is for research use only and is not approved for use in humans or in clinical diagnosis.   |
| RefSeq:                               | <a href="#">NP_000568</a>   |
| Locus ID:                             | 3557  |
| Cytogenetics:                         | 2q14.1  |
| Protein Families:                     | Druggable Genome, Secreted Protein  |



[View online »](#)