

## Product datasheet for **TB412511**

### Syntenin (SDCBP) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of SDCBP (NM_001007070), transcript variant 5, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC212511
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Syntenin (SDCBP) Mouse Monoclonal Antibody [Clone ID: OT12H6] (TA504796)
ACCN:	<a href="#">NM_001007070</a> , <a href="#">NP_001007071</a>
Synonyms:	MDA-9; MDA9; ST1; SYCL; TACIP18
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001007071</a>
Locus ID:	6386
Cytogenetics:	8q12.1
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)