

Product datasheet for **TB412674**

PTOP (ACD) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ACD (NM_001082487), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC212674
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PTOP (ACD) Mouse Monoclonal Antibody [Clone ID: OTI1A2] (TA504406)
ACCN:	<u>NM_001082487</u> , <u>NP_001075956</u>
Synonyms:	PIP1; PTOP; TINT1; TPP1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_001075956</u>
Locus ID:	65057
Cytogenetics:	16q22.1



[View online »](#)