

Product datasheet for **TB412712**

MAGI1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of MAGI1 (NM_001033057), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC212712
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	NM_001033057 , NP_001028229
Synonyms:	AIP-3; AIP3; BAIAP1; BAP-1; BAP1; MAGI-1; MAGI-1b; Magi1d; TNRC19; WWP3
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001028229
Locus ID:	9223
Cytogenetics:	3p14.1
Protein Pathways:	Tight junction



[View online »](#)