

## Product datasheet for **TB412837**

### **GABPB2 (GABPB1) CytoSection**

#### **Product data:**

<b>Product Type:</b>	CytoSections
<b>Description:</b>	Transient overexpression of GABPB1 (NM_181427), transcript variant gamma-3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	TrueORF Clone RC212837
<b>Tag:</b>	C-MYC/DDK
<b>Detection Antibodies:</b>	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
<b>Target Detection Antibodies:</b>	GABPB2 (GABPB1) Mouse Monoclonal Antibody [Clone ID: OTI7A7] (TA811268)
<b>ACCN:</b>	<a href="#">NM_181427</a> , <a href="#">NP_852092</a>
<b>Synonyms:</b>	BABPB2; E4TF1; E4TF1-47; E4TF1-53; E4TF1B; GABPB; GABPB-1; GABPB2; NRF2B1; NRF2B2
<b>Storage:</b>	Room Temperature, or 2-8°C for long term storage
<b>Stability:</b>	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
<b>Preparation:</b>	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
<b>Note:</b>	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
<b>RefSeq:</b>	<a href="#">NP_852092</a>
<b>Locus ID:</b>	2553
<b>Cytogenetics:</b>	15q21.2
<b>Protein Families:</b>	Transcription Factors



[View online »](#)