

## Product datasheet for **TB412964**

### DHLAG (CD74) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CD74 (NM_001025159), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC212964
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	DHLAG (CD74) Mouse Monoclonal Antibody [Clone ID: OTI1H3] (TA507339)
ACCN:	<u><a href="#">NM_001025159</a></u> , <u><a href="#">NP_001020330</a></u>
Synonyms:	DHLAG; HLADG; Ia-GAMMA; II; p33
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001020330</a></u>
Locus ID:	972
Cytogenetics:	5q33.1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Antigen processing and presentation



[View online »](#)