

Product datasheet for **TB413115**

FBXO22 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of FBXO22 (NM_012170), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC213115
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	FBXO22 Mouse Monoclonal Antibody [Clone ID: OTI1H1] (TA809766)
ACCN:	NM_012170 , NP_036302
Synonyms:	FBX22; FISTC1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_036302
Locus ID:	26263
Cytogenetics:	15q24.2
Protein Families:	Druggable Genome



[View online »](#)