

Product datasheet for **TB413159**

TAF2 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of TAF2 (NM_003184) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC213159
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	TAF2 Mouse Monoclonal Antibody [Clone ID: OTI2A5] (TA504588)
ACCN:	NM_003184 , NP_003175
Synonyms:	CIF150; MRT40; TAF2B; TAFII150
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_003175
Locus ID:	6873
Cytogenetics:	8q24.12
Protein Families:	Protease, Transcription Factors
Protein Pathways:	Basal transcription factors



[View online »](#)