

## Product datasheet for **TB413311**

### MSK2 / RSK-B (RPS6KA4) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RPS6KA4 (NM_001006944), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC213311
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_001006944</a></u> , <u><a href="#">NP_001006945</a></u>
Synonyms:	MSK2; RSK-B; S6K-alpha-4
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_001006945</a></u>
Locus ID:	8986
Cytogenetics:	11q13.1
Protein Families:	Druggable Genome, Protein Kinase, Transcription Factors
Protein Pathways:	MAPK signaling pathway, Neurotrophin signaling pathway



[View online »](#)