

Product datasheet for **TB413389**

ACVRL1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of ACVRL1 (NM_001077401), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC213389
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	ACVRL1 Mouse Monoclonal Antibody [Clone ID: OTI5C2] (TA809560)
ACCN:	NM_001077401 , NP_001070869
Synonyms:	ACVRLK1; ALK-1; ALK1; HHT; HHT2; ORW2; SKR3; TSR-I
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001070869
Locus ID:	94
Cytogenetics:	12q13.13
Protein Families:	Druggable Genome, Protein Kinase, Transmembrane
Protein Pathways:	Cytokine-cytokine receptor interaction, TGF-beta signaling pathway



[View online »](#)