

## Product datasheet for **TB413409**

### CTSA CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CTSA (NM_000308), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC213409
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u><a href="#">NM_000308</a></u> , <u><a href="#">NP_000299</a></u>
Synonyms:	GLB2; GSL; NGBE; PPCA; PPGB
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_000299</a></u>
Locus ID:	5476
Cytogenetics:	20q13.12
Protein Families:	Druggable Genome, Protease
Protein Pathways:	Renin-angiotensin system



[View online »](#)