

## Product datasheet for **TB413417**

### CD56 (NCAM1) CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NCAM1 (NM_001076682), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC213417
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	CD56 (NCAM1) Mouse Monoclonal Antibody [Clone ID: OT11G4] (TA506208)
ACCN:	<a href="#">NM_001076682</a> , <a href="#">NP_001070150</a>
Synonyms:	CD56; MSK39; NCAM
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001070150</a>
Locus ID:	4684
Cytogenetics:	11q23.2
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Cell adhesion molecules (CAMs), Prion diseases



[View online »](#)