

Product datasheet for **TB413863**

SAPAP3 (DLGAP3) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of DLGAP3 (NM_001080418) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC213863
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OT14C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
Target Detection Antibodies:	SAPAP3 (DLGAP3) Mouse Monoclonal Antibody [Clone ID: OT11D3] (TA811498)
ACCN:	NM_001080418 , NP_001073887
Synonyms:	DAP3; SAPAP3
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_001073887
Locus ID:	58512
Cytogenetics:	1p34.3



[View online »](#)