

Product datasheet for **TB414127**

Muscarinic Acetylcholine Receptor M4 (CHRM4) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CHRM4 (NM_000741) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC214127
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u>NM_000741</u> , <u>NP_000732</u>
Synonyms:	HM4; M4R
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_000732</u>
Locus ID:	1132
Cytogenetics:	11p11.2
Protein Families:	Druggable Genome
Protein Pathways:	Neuroactive ligand-receptor interaction, Regulation of actin cytoskeleton



[View online »](#)