

Product datasheet for **TB414238**

Caldesmon (CALD1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CALD1 (NM_033139), transcript variant 4, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC214238
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_033139</u> , <u>NP_149130</u>
Synonyms:	CDM; H-CAD; HCAD; L-CAD; LCAD; NAG22
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_149130</u>
Locus ID:	800
Cytogenetics:	7q33
Protein Pathways:	Vascular smooth muscle contraction



[View online »](#)