

Product datasheet for **TB414525**

CNG1 (CNGA1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of CNGA1 (NM_000087), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC214525
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_000087</u> , <u>NP_000078</u>
Synonyms:	CNCG; CNCG1; CNG-1; CNG1; RCNC1; RCNCa; RCNCalpha; RP49
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_000078</u>
Locus ID:	1259
Cytogenetics:	4p12
Protein Families:	Druggable Genome, Ion Channels: Cyclic nucleotide gated, Transmembrane



[View online »](#)