

## Product datasheet for **TB414570**

### BAIAP2 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of BAIAP2 (NM_006340), transcript variant 3, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC214570
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	BAIAP2 Mouse Monoclonal Antibody [Clone ID: OTI 1D9] (TA502981)
ACCN:	<u><a href="#">NM_006340</a></u> , <u><a href="#">NP_006331</a></u>
Synonyms:	BAP2; FLAF3; IRSP53; WAML
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u><a href="#">NP_006331</a></u>
Locus ID:	10458
Cytogenetics:	17q25.3
Protein Families:	Druggable Genome
Protein Pathways:	Adherens junction, Regulation of actin cytoskeleton



[View online »](#)