

## Product datasheet for **TB414617**

### RFC3 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RFC3 (NM_181558), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC214617
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	RFC3 Mouse Monoclonal Antibody [Clone ID: OTI3F6] (TA811874)
ACCN:	<a href="#">NM_181558</a> , <a href="#">NP_853536</a>
Synonyms:	RFC38
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_853536</a>
Locus ID:	5983
Cytogenetics:	13q13.2
Protein Families:	Stem cell - Pluripotency
Protein Pathways:	DNA replication, Mismatch repair, Nucleotide excision repair



[View online »](#)