

## Product datasheet for **TB414669**

### B3GALT5 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of B3GALT5 (NM_033170), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC214669
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<a href="#">NM_033170</a> , <a href="#">NP_149360</a>
Synonyms:	3-GalTase 5; B3GalT-V; B3GalTx; B3T5; beta-1; beta-3-Gx-T5; beta3Gal-T5; GLCT5
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_149360</a>
Locus ID:	10317
Cytogenetics:	21q22.2
Protein Families:	Transmembrane
Protein Pathways:	Glycosphingolipid biosynthesis - globo series, Glycosphingolipid biosynthesis - lacto and neolacto series, Metabolic pathways



[View online »](#)