

Product datasheet for **TB414670**

RAD17 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RAD17 (NM_133343), transcript variant 6, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC214670
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
ACCN:	<u>NM_133343</u> , <u>NP_579921</u>
Synonyms:	CCYC; HRAD17; R24L; RAD17SP; RAD24
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_579921</u>
Locus ID:	5884
Cytogenetics:	5q13.2
Protein Families:	Druggable Genome



[View online »](#)