

Product datasheet for **TB414699**

RASA1 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of RASA1 (NM_022650), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC214699
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	RASA1 Mouse Monoclonal Antibody [Clone ID: OTI4E2] (TA505935)
ACCN:	NM_022650 , NP_072179
Synonyms:	CM-AVM; CMAVM; CMAVM1; GAP; p120; p120GAP; p120RASGAP; PKWS; RASA; RASGAP
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_072179
Locus ID:	5921
Cytogenetics:	5q14.3
Protein Families:	Druggable Genome
Protein Pathways:	Axon guidance, MAPK signaling pathway



[View online »](#)