

## Product datasheet for **TB414701**

### IFIT3 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of IFIT3 (NM_001031683), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC214701
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	IFIT3 Mouse Monoclonal Antibody [Clone ID: OTI1G1] (TA500726)
ACCN:	<a href="#">NM_001031683</a> , <a href="#">NP_001026853</a>
Synonyms:	CIG-49; cig41; GARG-49; IFI60; IFIT4; IRG2; ISG60; P60; RIG-G
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_001026853</a>
Locus ID:	3437
Cytogenetics:	10q23.31



[View online »](#)