

Product datasheet for **TB415175**

PPP2R1B CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PPP2R1B (NM_002716), transcript variant 1, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215175
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PPP2R1B Mouse Monoclonal Antibody [Clone ID: OTI8B6] (TA811413)
ACCN:	<u>NM_002716</u> , <u>NP_002707</u>
Synonyms:	PP2A-Abeta; PR65B
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_002707</u>
Locus ID:	5519
Cytogenetics:	11q23.1
Protein Families:	Druggable Genome, Phosphatase, Transcription Factors
Protein Pathways:	Long-term depression, Oocyte meiosis, TGF-beta signaling pathway, Tight junction, Wnt signaling pathway



[View online »](#)