

Product datasheet for **TB415250**

ARFBP1 (HUWE1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of HUWE1 (NM_031407) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215250
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	<u>NM_031407</u> , <u>NP_113584</u>
Synonyms:	ARF-BP1; HECTH9; HSPC272; Ib772; LASU1; MRXST; MULE; URE-B1; UREB1
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<u>NP_113584</u>
Locus ID:	10075
Cytogenetics:	Xp11.22
Protein Families:	Druggable Genome
Protein Pathways:	Ubiquitin mediated proteolysis



[View online »](#)