

Product datasheet for **TB415266**

Kv4.2 (KCND2) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of KCND2 (NM_012281) in HEK293T cells paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215266
Tag:	C-MYC/DDK
Detection Antibodies:	Clone OTI4C5, Anti-DDK (FLAG) monoclonal antibody (TA50011-100)
ACCN:	NM_012281 , NP_036413
Synonyms:	KV4.2; RK5
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_036413
Locus ID:	3751
Cytogenetics:	7q31.31
Protein Families:	Druggable Genome, Ion Channels: Potassium, Transmembrane



[View online »](#)