

Product datasheet for **TB415331**

Kv beta 1 (KCNAB1) CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of KCNAB1 (NM_003471), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215331
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	Kv beta 1 (KCNAB1) Mouse Monoclonal Antibody [Clone ID: OTI7F12] (TA503971)
ACCN:	NM_003471 , NP_003462
Synonyms:	AKR6A3; hKvb3; hKvBeta3; KCNA1B; KV-BETA-1; Kvb1.3
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_003462
Locus ID:	7881
Cytogenetics:	3q25.31
Protein Families:	Druggable Genome, Ion Channels: Other



[View online »](#)