

## Product datasheet for **TB415348**

### PAX3 CytoSection

#### Product data:

Product Type:	CytoSections
Description:	Transient overexpression of PAX3 (NM_181460), transcript variant PAX3H, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215348
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	PAX3 Mouse Monoclonal Antibody [Clone ID: OTI4D1] (TA809630)
ACCN:	<a href="#">NM_181460</a> , <a href="#">NP_852125</a>
Synonyms:	CDHS; HUP2; WS1; WS3
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	<a href="#">NP_852125</a>
Locus ID:	5077
Cytogenetics:	2q36.1
Protein Families:	Adult stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Transcription Factors



[View online »](#)