

Product datasheet for **TB415425**

NR2E3 CytoSection

Product data:

Product Type:	CytoSections
Description:	Transient overexpression of NR2E3 (NM_014249), transcript variant 2, in HEK293T cells, paraffin embedded controls for ICC/IHC staining
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	TrueORF Clone RC215425
Tag:	C-MYC/DDK
Detection Antibodies:	DDK Rabbit monoclonal antibody, recognizing both N- and C-terminal tags (TA592569)
Target Detection Antibodies:	NR2E3 Mouse Monoclonal Antibody [Clone ID: OTI3D4] (TA806203)
ACCN:	NM_014249 , NP_055064
Synonyms:	ESCS; PNR; rd7; RNR; RP37
Storage:	Room Temperature, or 2-8°C for long term storage
Stability:	Blocks are guaranteed for a year from the date of receipt if proper storage instructions were followed.
Preparation:	HEK293T cells were transiently transfected with TrueORF cDNA plasmid. Transfected cells were cultured for 48hrs. After harvesting, the cultured cells were fixed in formalin & dehydrated before embedding in paraffin.
Note:	This product is for research use only and is not approved for use in humans or in clinical diagnosis.
RefSeq:	NP_055064
Locus ID:	10002
Cytogenetics:	15q23
Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors



[View online »](#)